User's Manual
ENGLISH

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USER MANUAL

DESCRIPTION

 is a non-absorbable, permanent-action, tissue bulking substance, for the endoscopic treatment of vesicoureteral reflux (VUR). It is easily and manually injected.

 consists of particles of polyacrylate polyalcohol immersed in a glycerol and saline solution carrier. It has a very high molecular mass (~10 million Daltons) and it comes in the form of sterile, pyrogen-free particles that are highly deformable by compression, and that may be injected with 23-gauge needles.

Once implanted, no local, regional or distance migration has been observed. The carrier is a 40% glycerol solution. Once implanted, it is eliminated by the reticuloendothelial system without metabolizing and excreted through the kidneys, while the particles remain to ensure permanent bulking.

INSTRUCTIONS FOR USE AND MODE OF ACTION

The physicochemical properties of  make it suitable for the endoscopic treatment of vesicoureteral reflux (VUR).

When  is implanted in the ureterovesical junction, the material enlarges the volume of the area and corrects the anatomy of the meatus and the distal ureter, preventing urine to return to the ureter after having been stored in the bladder.

HOW SUPPLIED

 (code BAR- 1J) is supplied as follows:

- 1 syringe (Code: BARI-1J) with stopper, containing 1 ml of material.

 is provided sterile and ready-to-use. The syringe is sterilized with moist heat.

VANTRIS VUR can be injected using either a flexible injection needle or a metal semi-rigid injection needle.

It is recommended to use the injection needles developed by Promedon, with the following features:

- 1 flexible injection needle (Code 50 F), 5 Fr x 23 G x 305 mm.

Or alternatively:

- 1 metal semi-rigid injection needle (Code RIN), 3.6 Fr x 22 G x 350 mm.

HANDLING AND STORAGE PRECAUTIONS

-  is provided sterile and pyrogen-free.
- : The package consists of a cardboard box containing the syringe.

The syringe is double-wrapped and protected by an aluminum pouch which should be removed from the surgical field. IF ANY OF THE 3 POUCHES IS DAMAGED, DO NOT IMPLANT .

Packaging

 is supplied in a box containing the syringe in a triple pouch; the external one is made of aluminum.

HANDLING AND STORAGE PRECAUTIONS

 have been designed to be used only ONCE. Therefore, DO NOT REUSE or RSTERILIZE.

IF ANY POUCH CONTAINING VANTRIS VUR COMPONENTS WERE DAMAGED, DO NOT IMPLANT.

IMPORTANT: Only the internal pouch can enter the sterile area.

Storage

 should be stored under the following conditions:

- TEMPERATURE: Room temperature.

DO NOT USE AFTER THE EXPIRATION DATE ON THE PACKAGING.

OPERATING ROOM CONDITIONS

In general, all operating rooms comply with the necessary conditions to perform an implant surgery. However, the following parameters should be specially verified:

- Operating room asepsis.
- Adequate training of the operating room staff

PATIENT INFORMATION

Surgeons are responsible for informing patients or their representatives on

the possible complications associated to the VANTRIS VUR implant surgery, in advance.

CONTRAINDICATIONS

- Ureterocele
- Voiding dysfunction
- Kidney dysfunction
- Paraureteral (Hutch) diverticulum
- Urinary tract infections
- Uterovesical junction obstruction

PRECAUTIONS

The technique and the instruments necessary for injecting  should follow the same precautions as in other similar urological interventions, endoscopy practices and endourological procedures.

 should be used by qualified surgeons with experience in endoscopy practices and trained in subureteral injection procedures.

Antibiotic prophylaxis following the treatment with  may be advisable and its application depends on the surgeon's criteria.

WARNINGS

Tissues with exacerbated scarring processes and/or with fibrosis can compromise the ideal bulking.

 should not be injected intravascularly since it may cause vascular occlusion.

COMPLICATIONS

The possible complications associated with the use of  should be discussed with the patient's parents prior to the surgery.

Complications associated with the implantation of an injectable agent into the vesicoureteral junction include those common to this types of procedures:

- Temporary post-injection dilatation of the upper urinary tract.
- Obstruction of the vesicoureteral junction.
- Dysuria
- Hematuria
- Infection
- Technique failure

Promedon requires that all surgeons inform the Company of any complication observed with the use of .

IMPLANT PROCEDURE

In order to perform the procedure, the following components are necessary:

• 1 syringe (Code: BARI-1J): The quantity to be used will depend on each particular case.

• 1 flexible injection needle or alternatively 1 metal semi-rigid injection needle. It is recommended to used the injection needles developed by Promedon, with the following features:

- 1 flexible injection needle (Code 50 F), 5 Fr x 23 G x 305 mm.
Or alternatively:
- 1 metal semi-rigid injection needle (Code RIN), 3.6 Fr x 22 G x 350 mm.

• 1 cystoscope with a working channel of 4 Fr or more, according to the needle chosen.

1.  is injected under general anesthesia in children.
2. The free flow through the injection needle is verified using saline solution.
3. The syringe is connected to the injection needle passing material through it until the substance appears at the needle tip. Insert the injection needle into the cystoscope.
4. Prior to the injection, the meatus should be observed with different volumes of bladder filling in order to choose the ideal situation.
5. The submucosa of the bladder is punctured at the six o'clock position, 3 mm under the ureteral meatus, 4 to 5 mm depth.
6. The material is then slowly injected until the ureteral wall is adequately modified (bulkiness). More punctures may be performed until the desired effect is reached.
7. Once the injection is completed, the needle is kept in its position for 30 seconds.
Then, the needle is removed.
8. Once the procedure has been completed, the bladder is emptied, the optic removed, and an optional voiding cystoureterogram (VCUG) may be performed to verify reflux repair.

NOTES

Consider that 0.4 ml of the product remains in the injection needle, which can be used by pushing it down with the same syringe with sterile physiological solution or sterile distilled water.
A ureteral catheter may be necessary where the meatus presents difficulties for injection.

Re-implantation needs: If after performing a voiding cystoureterogram (VCUG) three months after the treatment the persistence of reflux is proven, Vantris VUR can be re-injected until it is eliminated.

SYMBOLS USED IN THE LABELS

CATALOGUE NUMBER

BATCH NUMBER

CAUTION

EXPIRATION DATE

DO NOT REUSE

EXPIRATION DATE

STERILE. STERILIZATION METHOD:
MOIST HEAT

MANUFACTURER

AUTHORIZED REPRESENTATIVE
IN THE EUROPEAN COMMUNITY

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